

Form B03

Scheduled Maintenance Work Order

Quantum Medical Solutions sdn bhd
transcending boundaries, transforming life

Format Ref. - QMS/TSD-003 Rev. 00

Work Order No	PW0371951	Scheduled Month	FEBRUARY 2018
Work Order Date	1/2/18	Completed Date	15/2/18
Clinic Name	KD BATU MANIKAR	Clinic Code	WPL007
BE No	WPL000504	District	W.P. LABUAN
BE Category	FETAL, HEART DETECTORS	WO Assigned to	che Muhammed
Ownership	<input checked="" type="checkbox"/> Existing Equipment	<input type="checkbox"/> New & Purchased Equipment	
BE Condition	<input type="checkbox"/> Under Warranty	<input type="checkbox"/> BER Proposed	
Work Order Type	<input checked="" type="checkbox"/> Preventive Maintenance (PM) <input type="checkbox"/> Routine Inspection (RI)	<input type="checkbox"/> Third Party Calibration (TPC) <input type="checkbox"/> Statutory Certification (SC)	

BE Third Party Calibration / Statutory Certification Details

Company Name	NA	Cal / Cert Date	NA
Contact Number		Cal / Cert Expiry Date	

Action Taken

↳ PPM was carried out as per check.

Schedule Maintenance Execution Details

SI No	QMS Engineer / Technician Name	Date	Start Time	End Time
	MAZIANA BT YUNUS	15/2/18	1015	1030

Customer Remarks

--

Engineer / Technician Signature Name: <i>MAZIANA BT YUNUS</i> Date: 15/2/18	Customer Signature Name: <i>MAZIANA BT YUNUS</i> Designation: JURURAWAT W. SARAWAK Date: 15.02.2018 Stamp:
---	--


For Internal Use only


First Verification
QMS Circle Incharge

JULIUS LIANSUN
BIOMEDICAL ENGINEER, QMS
019-3620179

Final Verification
QMS State Incharge

DICKY LEE
SABAH STATE MANAGER
QUANTUM MEDICAL SOLUTIONS SDN BHD

 Quantum Medical Solutions sdn bhd <small>Powering homecare, transforming life</small>	Quantum Medical Solutions Sdn Bhd BEMS Planned Preventive Maintenance Checklist Fetal Heart Detectors, Ultrasonic BE CODE : 11-696	CHECKLIST NO: 070 REV.000					
PART 1 ASSET DETAILS							
WORK ORDER NO ▶ pw 0371951		ASSET NO ▶ WFL000504					
MANUFACTURER ▶ Sonicare one		MODEL ▶ N1A					
FREQUENCY ▶ 3 MONTHLY () 6 MONTHLY () 12 MONTHLY (✓)		PPM HOURS ▶ 0.25					
PART 2 SPECIAL PRECAUTION							
If there is evidence of body fluid contamination, submit the device for cleaning and decontamination before inspecting it. Wear appropriate Personnel Protection Equipment (PPE) during work. Wear grounded electrostatic wristband when handling PCB or electronic components. Refer to the safety procedure for additional precautions and guidance as per manufacturer guidelines. Make sure the test equipment used are duly calibrated.							
PART 3 TEST APPARATUS							
<i>Tick (✓) where appropriate</i>							
NO	ASSET NO	DESCRIPTION	SERIAL NO	CALIBRATION DUE ON			
	N/A	ELECTRICAL SAFETY TESTER	NA	NA			
PART 4 QUALITATIVE TASKS							
<i>Tick (✓) where appropriate</i>							
	PASS	FAIL	NA		PASS	FAIL	NA
1 Chassis - verify physical integrity, cleanliness and condition	(✓)	()	()	7 Circuit Breaker/ Fuse - verify integrity of external circuit breaker and/or rating of external fuse	(✓)	()	()
2 Mount/ Fasteners - verify physical integrity	(✓)	()	()	8 Fittings/ Connectors - check all fittings/connectors	(✓)	()	()
3 Casters/Brakes - if mounted, verify physical integrity	()	()	(✓)	9 Controls/ Switches/ Keypad - verify proper operation of controls	(✓)	()	()
4 Power Cord - verify proper insulation and integrity	()	()	(✓)	10 Indicators/ Displays - verify proper illumination and operation	(✓)	()	()
5 Strain Relief - verify physical integrity at both ends of line cord	()	()	(✓)	11 Alarm/ Audible Signal - verify operation	(✓)	()	()
6 Transducers/Cables - verify integrity and condition	()	()	(✓)				
PART 5 PREVENTIVE MAINTENANCE TASKS							
<i>Tick (✓) where appropriate</i>							
	DONE	NOT DONE **	NA				
1 Cleaning the inside with a vacuum cleaner	()	()	(✓)	Notes: * For all Parts, NA is defined as NOT APPLICABLE ** If you have ticked 'NOT DONE', then justify in Part 8 *** Choose Whichever Applicable			
2 Clean Transducers/Cables	(✓)	()	()				
3 Check/Replace battery ***	(✓)	()	()				



REMS Planned Preventive Maintenance Checklist

Fetal Heart Detectos,Ultrasonic

BE CODE : 11-696

CHECKLIST NO: CL-070
REV.000

WORK ORDER NO ▶ JWO-371951

PART 6 QUANTITATIVE TASKS

Tick (\checkmark) where appropriate

[illegible]

PART 7 ELECTRICAL SAFETY TEST

ELECTRICAL SAFETY TEST, (attach report)


(In accordance to IEC 60601)

☐ PASS☐ FAIL☒ NA

PART 8 NOTES

MAZIANA BT YENUS
JURURAN
No. LIM : 9691

15.02.18



☐ CORRECTIVE MAINTENANCE REQUIRED

☒ FUNCTIONING

☐ NOT FUNCTIONING

WORK ORDER NO. NA

NEXT PPM DATE ► 14/2/14

PPM has been performed in accordance to the checklist and the equipment is functioning to the intended purpose.

COMPLETED BY:

DATE :

15/2/18